

Durasul® Highly Crosslinked Polyethylene
Advanced Tribology for Total Knee and Hip
Arthroplasty that Resists Wear and Aging



Durasul Highly Crosslinked Polyethylene

Advanced Tribology for Total Knee and Hip Arthroplasty



Highly crosslinked polyethylene acetabular components enhance the surgeon's ability to provide patients with a long-term solution.



Durasul Natural-Knee[®] tibial and patellar components were tested extensively under both normal and adverse conditions.

A breakthrough in total joint replacement, Durasul highly crosslinked polyethylene directly addresses the problem of wear that can lead to osteolysis and implant loosening.

Durasul's properties are stable over time, as demonstrated by testing after accelerated aging. Thus, Durasul has the potential to extend the implant's durability, providing a long-term solution for joint replacement.

Why compromise?

The three key requirements for polyethylene in joint replacement include:

- resistance to wear,
- resistance to aging, and
- device-suitable mechanical properties.¹

Durasul was designed and tested to fulfill all three.

A Leader in Innovation

1998 In collaboration with Dr. William Harris at Massachusetts General Hospital (MGH) and MIT, we developed Durasul highly crosslinked polyethylene.

First Durasul acetabular component implanted in Sweden.

1999 Durasul is the 1st highly crosslinked acetabular system available in U.S. market.

Market surveillance studies begin at over 300 centers worldwide.

2000 Durasul enables introduction of the 1st Large Diameter Head System.

2001 Durasul tibial inserts introduced to Natural-Knee[®] System - 1st knee system to use highly crosslinked polyethylene.

2002 Another 1st: Durasul patella introduced.


The results of in-vitro testing have not been shown to correlate with clinical wear mechanisms.

¹ Muratoglu OK, Bragdon CR, O'Connor DO, Jasty M, Harris WH. A novel method of cross-linking ultra-high-molecular-weight polyethylene to improve wear, reduce oxidation, and retain mechanical properties. *Journal of Arthroplasty*. 2001; 16(2):149-60

Durasul Highly Crosslinked Polyethylene


Resistant to Wear and Aging without Compromising Mechanical Properties

Wear Rate vs. Radiation Dose



Bi-directional pin-on-disk wear test. Durasul is irradiated at 9.5 Mrad to maximize wear benefits⁵ and maintain mechanical properties.

Aging Effect on Mechanical Properties^{8,9}



Unlike conventional polyethylene, Durasul meets ASTM and ISO minimums in both the aged and unaged conditions.

Durasul is highly crosslinked to resist wear.

Crosslinking is the formation of chemical bonds between adjacent molecular chains, making them more difficult to separate and improving wear properties. This has been demonstrated clinically^{2,3}, as well as through simulator testing.^{1,4}

Free radicals removed to prevent aging.


During irradiation, some free radicals form crosslinks while others do not. The remaining free radicals may react with oxygen over time, leading to embrittlement, reduction in mechanical properties, and increased wear.^{6,7} This change in properties over time is an example of “aging.”

To remove the remaining free radicals, Durasul is heated above its melt temperature following crosslinking, allowing residual free radicals to recombine. This creates additional crosslinks while ensuring the stability of Durasul over time.

Device testing ensures appropriate mechanical properties for Durasul components.

Centerpulse Orthopedics, in conjunction with MGH, pioneered rigorous device-specific biomechanical and simulator tests for Durasul tibial, patellar, and acetabular components. In all cases, the Durasul components met or surpassed the performance of conventional polyethylene components.

A unique manufacturing process designed from the latest in materials technology




1. Oonishi H. Long term clinical results of THR: clinical results of THR of an alumina head with cross-linked UHMWPE cup. *Orthop Surg Traumatol.* 1995; 38:1225.
2. Wroblewski BM, Siney PD, Dowson D, Colling SN. Prospective clinical and joint simulator studies of a new total hip arthroplasty using alumina ceramic heads and cross-linked polyethylene cups. *Journal of Bone and Joint Surgery (Br).* 1996; 78:280.
3. McKellop H, Shen F, Lu B, Campbell P, Salovey R. Development of an extremely wear-resistant ultra high molecular weight polyethylene for total hip replacements. *Journal of Orthopaedic Research.* 1999; 17:157-67.
4. Muratoglu OK, Bragdon CR, O'Connor DO, Jasty M, Harris WH, Gul R, McGarry F. Unified wear model for highly crosslinked ultra-high molecular weight polyethylenes (UHMWPE). *Biomaterials.* 1999; 20(16):1463-70.
5. Fisher J, Reeves EA, Isaac GH, Saum KA, Sanford WM. Comparison of the wear of aged and non-aged ultrahigh molecular weight polyethylene sterilized by gamma irradiation and by gas plasma. *Journal of Material Science: Materials in Medicine.* 1997; 8(6): 375-8.
6. Bohl JR, Bohl WR, Postak PD, Greenwald AS. The Coventry Award. The effects of shelf life on clinical outcome for gamma sterilized polyethylene tibial components. *Clin Orthop.* 1999; Oct(367):28.
7. Muratoglu OK, Bragdon CR, O'Connor D, Perinchieff RS, Estok DM, Jasty M, Harris WH. Larger Diameter Femoral Heads Used in Conjunction with a Highly Cross-linked Ultra-High Molecular Weight Polyethylene. *J Arthroplasty.* 2001; 16(8) Suppl 1:24-30.
8. DiMaio WG, Saum KA, Lily WB, Moore WC. Effect of Sterilization Type and Accelerated Aging on Physical Properties of UHMWPE. *24th Annual Meeting of the Society for Biomaterials.* San Diego, April 22-26, 1998.

Durasul for the Knee

Addressing Delamination and Wear

These simulator studies have been correlated to implant retrievals to substantiate Durasul's superior wear performance over conventional polyethylene.^{11,12}



We pioneered new standards for tibio-femoral and patello-femoral testing where none previously existed.

Patello-Femoral Simulator Results

	Retrieved conventional polyethylene, after 4.5 years <i>in vivo</i> . Note severe delamination on the lateral edge.
	Aged conventional polyethylene, after 1 million cycles on simulator. Damage mode is similar to that in retrievals.
	Aged Durasul, after 2 million cycles on simulator. No delamination observed.

¹⁰ Muratoglu OK, O'Connor DO, Bragdon CR, Jasty M, Perinchief RS, Savadove T, Harris WH, Krevolin J, Shen M. Wear Evaluation of Tibial Inserts made from Conventional and Durasul UHMWPE through simulated Gait. *White paper #1001-01-042*. 2001.

¹¹ Muratoglu OK, Burrough BR, O'Connor DO, Bragdon CR, Perinchief RS, Sargent MR, Krevolin J, Shen M, Estok II DM, Freiburg AA, Rubash HE, Jasty M, Harris WH. Comparison of Conventional and Durasul UHMWPE on a new In Vitro Patella simulator. *White paper #1001-01-048*. 2001.

¹² Krevolin J, Shen M, Bailey A, Prybyla R, Muratoglu OK, Burroughs BR, Delozio K, Wyss U, Harris WH. Comparison of Conventional Versus Durasul UHMWPE for the Patella: A New In Vitro Simulation of Stair Climbing. *White paper #1001-01-050*. 2001.


Durasul for the Hip

Larger Heads for Increased Stability

The unique combination of Durasul and large heads addresses the two most common complications following hip replacement: wear-induced osteolysis and postoperative dislocation of the hip.^{13,14}

Durasul liners for use with the Converge® Acetabular System demonstrated immeasurable wear in hip simulator tests after 27 million cycles.¹⁵ In fact, the wear of Durasul liners is independent of head size, while wear of conventional polyethylene liners increases dramatically as head size increases.¹⁵ Thus, femoral heads with diameters up to 44 mm may be paired with Converge Durasul liners. Larger diameter heads more closely replicate the patient's anatomy, enhance stability, and increase ROM.

Durasul acetabular liners performed well in biomechanical tests, including: fatigue¹⁶, micromotion¹⁷, push-out/lever-out¹⁸, and contact stress analysis¹⁹ with inserts down to 3 mm thick.



Converge hemispherical acetabular cup with Durasul liner and 44-mm CoCr head on a Natural-Hip™ stem.



Tested under normal and adverse conditions

Durasul inserts are tested within their shells ensuring anatomical testing of specific design features.

13 Turner RS. Postoperative total hip prosthetic femoral head dislocations. Incidence, etiologic factors, and management. *Clinical Orthopedics*. 1994; 301:196-204.

14 McCollum DE, Gray WJ. Dislocation after total hip arthroplasty. Causes and prevention. *Clinical Orthopedics*. 1990; 261:159-70.

15 Muratoglu OK, Bragdon CR, O'Connor D, Perinchieff RS, Estok DM, Jasty M, Harris WH. Larger Diameter Femoral Heads Used in Conjunction with a Highly Cross-linked Ultra-High Molecular Weight Polyethylene. *J Arthroplasty*. 2001; 16(8), Suppl 1:24-30.

16 Bailey A. Application of 3mm Durasul Polyethylene with the Inter-Op Acetabular System: Evaluation of Fatigue Performance Characteristics of the Size 41 x 28 mm Acetabular System Under Cyclic Physiological Loads and Torque. *Centerpulse Orthopedics Technical Report #99028G*. 2001.

17 Bailey A. Application of 3mm Durasul Polyethylene with the Inter-Op Acetabular System: Evaluation of Micromotion Performance Characteristics of the Size 41 x 28 mm Acetabular System Under Cyclic Physiological Loads and Torque. *Centerpulse Orthopedics Technical Report #99029G*. 2000.

18 Kralovic BJ. Evaluation of Shell-Insert Push-In, Push-Out, Lever-Out and Torque-Out of Converge Cluster-Hole 53 mm Shells and 28 x 54 mm Durasul Inserts. *Centerpulse Orthopedics Technical Report #01034G*. 2001.

19 Kurtz S, Giddings V. Stress analysis of the Inter-Op Durasul System. *Exponent Inc*. January 4, 2000.

Zimmer

9900 Spectrum Drive
Austin, Texas 78717

512.432.9900
877.768.7349
fax: 512.432.9014



zimmer
Confidence in your hands™

1001-00-013 Rev. 10/2003
© 2003 Zimmer
All rights reserved.

Knees

- Apollo® Knee System
 - Classic condylar knee replacement system
- Durasul® Tribological System
 - Highly crosslinked polyethylene that resists wear and aging
- Natural-Knee® System
 - Anatomic design for superior clinical results
- UniSpacer™ Knee System
 - No bone cuts. No compromises.

Severe Revision/Limb Salvage

- MOST Options™ System
 - Modular knee and hip options for severe bone loss, trauma and revision

Hips

- Alloclassic® (Zweymüller™) Hip
 - Classic proven design with superior clinical results
- Allofit™ Acetabular Cup System
 - Unique Ridgelock™ surface designed for easy implantation and stability
- Apollo® Hip System
 - Designed for optimal results with low-demand patients
- APR® Anatomical Hip System
 - The anatomic solution for bone matching
- CLS™ (Spotorno™) Hip System
 - The standard of proximal press-fit design
- Converge® CST™ Porous Acetabular Cup System
 - Where technology and experience meet
- Durasul® Tribological System
 - Highly crosslinked polyethylene that resists wear and aging
- FracSure™ Hip System
 - A classic design for hip fractures
- Metasul® Metal-on-Metal Acetabular System
 - Over 15 years of clinical results & 200,000 implantations worldwide
- MS-30™ Hip
 - A highly polished cemented stem
- Natural-Hip™ System
 - A comprehensive system with a natural approach
- Precedent™ Revision Hip System
 - A better solution for revision hips
- SL Revision™ Hip System
 - A stable revision design with extensive sizes

Upper Extremities

- Anatomical™ Shoulder System
 - Multiple adjustments of inclination & retroversion with the potential for precisely restored anatomy
- GSB® Elbow System
 - A nonconstrained design with 21 years of clinical results
- Select® Shoulder System
 - TSA and fracture management with offset head options