



**Drill Guide**  
**3.5mm Screw**  
**00-2345-035-25**

**Assembly/Disassembly**  
**Instructions**



## Disassembly Instructions for Drill Guide 3.5mm Screw

1. Pull the Guide Assembly **2** firmly to disengage it from the Handle Assembly **1**.
2. Pull the Guide Assembly **3** firmly to disengage it from the Handle Assembly **1**.
3. End of disassembly.

## Assembly Instructions for Drill Guide 3.5mm Screw

1. Push the Guide Assembly **2** firmly inside the Handle Assembly **1** until the spherical ball of the Guide Assembly **2** engages and free to rotate in the groove of the Handle Assembly **1**.
2. Push the Guide Assembly **3** firmly inside the Handle Assembly **1** until the spherical ball of the Guide Assembly **3** engages and free to rotate in groove of the Handle Assembly **1**.
3. End of assembly.



+H124970000506001/\$100610F100

97-0000-506-00 1003-T27 Printed in USA ©2010 Zimmer, Inc.



**zimmer** | trauma