




**Percutaneous Strut  
Cannula (1.6mm)  
00-2360-025-16**

**Assembly/Disassembly  
Instructions**



**2** Cannula Nut

Cannula Shaft **1**

## Disassembly Instructions for Percutaneous Strut Cannula

1. Unthread Cannula Nut **2** from Cannula Shaft **1**.
2. End of disassembly.

## Assembly Instructions for Percutaneous Strut Cannula

1. Thread Cannula Nut **2** to Cannula Shaft **1**.
2. End of assembly.



+H124970000508001/\$100610F10Q

97-0000-508-00 1003-T25 Printed in USA ©2010 Zimmer, Inc.